


IP-1020

Tabletop pulse oximeter

Professional medical care, real time operation system



2-Parameter vital signs
(SpO2+NIBP)



Compact & light weight



High sensitive NIBP technology
(neonate to adult)



Advanced
SpO2 technology





72 hours trend



Perfusion range:
0.05%-20%

Offer various mode

Oxygen saturation, pulse rate measurement (optional high sensitive NIBP)



Plethysmographic mode



SpO2 averaging



Trend data



Special motion tolerance algorithm

Measure the exact value of pure oxygen saturation by removing noise.

| Manufacturer



Infunix Technology Co., LTD

T +82-2-3448-2800

F +82-2-546-8346

E sales@infunix.com

| Sales



GENORAY Co.,Ltd.

T +82-31-5178-5500

E genoray@genoray.com

W www.genoray.com

GENORAY America Inc.

T +1-855-436-6729

E inquiry@genorayamerica.com

W www.genorayamerica.com

GENORAY EU GmbH

T +49-30-50969498

E info@genorayeuropa.com

W www.europe.genoray.com

GENORAY Japan

T +81-45-620-4971

E Info@genorayjapan.co.jp

W www.genorayjapan.co.jp

GENORAY Türkiye

T +90-216-463-57-37

+90-216-463-57-53

E genorayturkiye@genoray.com

GENORAY Shanghai

T +86-021-6785-8538

W www.genoray.com.cn